

FUSIONFORM DESKTOP™

OVERVIEW

FusionForm Desktop is an innovative, new product that captures and transforms handwritten data into digital data that is easily integrated into EMR (Electronic Medical Records) systems. Healthcare professionals can capture patient information in a natural and familiar manner and still satisfy the healthcare industry's need for digital data.

Using a digital pen, desktop software and their own medical forms printed on digital paper, healthcare professionals can instantly turn their free-form clinical notes, diagnoses, prescriptions, and test requisitions into digital data. With a simple interface for the user to instantly view and verify the digital results of handwriting recognition, FusionForm Desktop provides the accessibility of pen and paper and the verified data integrity that high quality healthcare requires.



TECHNICAL SPECIFICATIONS

System Requirements

Microsoft® Windows® 2000 or Microsoft Windows XP

Equipment Certified/Preferred

Anoto-enabled Logitech io2™ Digital Pen or Hitachi Maxell® Penit
FusionForm Desktop requires USB support

Export Formats Supported

XML, JPG and RTF

Bluetooth Technical Info (Bluetooth stack versions supported)

Widcomm™ Bluetooth stack version 1.2 or higher

Wireless Technologies

FusionForm has been designed to take advantage of the latest technologies and protocols including Bluetooth and Wi-Fi®.

Battery Life Estimates/Charging Digital Pens/User Experience

Digital pens: Under normal working conditions most digital pens hold their charge quite well. For best results, recharge when not in use in order to maintain charge.

Other Technical Information

Digital Pens: Replacement ballpoint refills for digital pens are available directly from the manufacturers or their resellers.

KEY FEATURES

- Departmental forms automation made easy
- Multiple healthcare users can contribute to a single form
- Desktop software can be easily integrated with healthcare EMRs and commercial databases
- Electronic forms images can be shared with others as “view only”
- Flexibility to incorporate existing workflow processes
- FusionForm allows users to select the level of forms automation they need based on the stage of technology deployment they have
- Unique dot pattern “fingerprints” improve patient safety, reduce liability, enhance data integrity and eliminate duplicate records

ABOUT SATORI LABS

Founded in 2003, Satori Labs, Inc. is a leader in forms automation solutions. Satori Labs fuses digital pen and paper, mobile and wireless technologies with Satori Labs' high-performance forms automation software to create the next generation of connected information management products for the mobile professional.

More information on Satori Labs technology is available at: www.satorilabs.com.