

FUSIONFORM MOBILE™

OVERVIEW

FusionForm Mobile is an innovative, new product that recognizes and captures handwritten information from customized medical forms and instantly transfers it into a PDA as digital data. Once in the FusionForm enabled PDA, this information is easily transmitted into an Electronic Medical Record (EMR) or other back-end database archive. Home healthcare workers and other field professionals can seamlessly enter data into a digital system using a familiar pen and paper format—without dictating, transcribing or key-boarding.

By using FusionForm Mobile, a digital pen, and their own forms printed on digital paper, mobile clinicians can easily and automatically convert free-form clinical notes, sketches, diagnoses, prescriptions or any other handwritten material into electronic data. A simple interface allows users to verify data at the point of entry on a PDA, and then either store it or wirelessly transmit it to a back-end.



TECHNICAL SPECIFICATIONS

System Requirements

Microsoft® Windows® Mobile 2003 SE operating system

Equipment Certified/Preferred

Hewlett-Packard® iPAQ PocketPC™ running Windows Mobile 2003 SE operating system
 Anoto-enabled Nokia SU-1B Digital Pen, Logitech io2™ Digital Pen, or Hitachi Maxell® Penit

Export Formats Supported

XML, JPG and RTF

Bluetooth Technical Info (Bluetooth stack versions supported)

Widcomm™ Bluetooth stack version 1.2 or higher

Wireless Technologies

FusionForm has been designed to take advantage of the latest technologies and protocols including Bluetooth and Wi-Fi®.

Battery Life Estimates/Charging Digital Pens/User Experience

PDAs: Please refer to manufacturer's product documentation for specific battery life guidelines.

Digital pens: Under normal working conditions most digital pens hold their charge quite well. For best results, recharge when not in use in order to maintain charge.

Other Technical Information

PDAs: Please refer to manufacturer's documentation for specific guidelines for proper maintenance and usage.

Digital Pens: Replacement ballpoint refills for digital pens are available directly from the manufacturers or their resellers.

KEY FEATURES

- Offers users true mobility when capturing client information
- FusionForm makes departmental forms automation easy
- FusionForm mobile can be easily integrated with healthcare EMRs and commercial databases
- Flexibility to incorporate existing workflow processes

ABOUT SATORI LABS

Founded in 2003, Satori Labs, Inc. is a leader in forms automation solutions. Satori Labs fuses digital pen and paper, mobile and wireless technologies with Satori Labs' high-performance forms automation software to create the next generation of connected information management products for the mobile professional.

More information on Satori Labs technology is available at: www.satorilabs.com.